


# BESIDE CABINETS

## Mesas de Noche

HF-01  
White-Cream



Main Material: Cold reduced metal sheet  
(0.03937" thickness)  
Metal color: White  
ABS Tabletop  
ABS Drawer and front door  
Extendable table  
Internal ABS shelf  
Size: 19.5"x20"x28"  
Castor size: 2"

HF-03  
Wooden Trim



Main Material: Cold reduced metal sheet  
(0.03937" thickness)  
Metal color: Beige  
ABS Tabletop  
ABS Drawer and front door  
Extendable table  
Internal ABS shelf  
Size: 19.5"x20"x28"  
Castor size: 2"